

RECORDS CHECKLIST

Single Dentate Arch Surgery: Opposing Edentulous

- Photographs** - Full face full smile photo, plus profile photo if Class II or III. Retracted photographs front, left side, right side IN occlusion, see below. We suggest capturing photos before taking impressions to keep the lips and teeth clean.
- Dentate surgery arch** - Master cast—must capture ALL land areas. Must capture full palate if spare denture is to be ordered.
- Edentulous opposing arch** - Impression capturing the teeth of the opposing dentition—a crown and bridge impression of the denture.
- Bite registration** - Suggestion: 3 small bites to tripod the bite registration rather than full arch bite
- CBCT** - Scan patient with the denture seated in occlusion. Export and upload multi-file .dcm. Zip the folder - should contain 200-400 individual .dcm files.

Photographs

Full Face / Full Smile + left / right / straight retracted in occlusion

