

## RECORDS CHECKLIST

### Single Edentulous Surgery: Opposing Dentate

- Photographs** - Full face full smile photo, plus profile photo if Class II or III. Retracted photographs front, left side, right side IN occlusion, see below. We suggest capturing photos before taking impressions to keep the lips and teeth clean.
- Edentulous surgery arch** - Use existing denture as a dual scan appliance. Denture must fit perfectly to the tissue and exhibit ideal tooth position. NO soft relines. Must use either hard reline if fit is not ideal, or reline with Blue Mousse impression material for both scans. Place 6 scan markers on the denture ([www.suremark.com](http://www.suremark.com)) or place 6 random gutta percha markers in the pink area of the labial flange.
- Dentate opposing arch** - Master cast must capture ALL land areas. Must capture full palate if spare denture is to be ordered.
- Bite registration** - Suggestion: 3 small bites to tripod the bite registration rather than full arch bite.
- CBCT** - Dual Scan Technique
  - \* *Dual Scan definition: one scan of the denture seated in the mouth IN occlusion. One scan of the denture alone sitting on foam or the chin rest of the scanner, not on plastic. Do not scan the patient alone w/o the denture.*
  - \* *Two scans are sent: 1) scan of maxillary arch 2) scan of mandibular arch. Both scan must be captured out of occlusion*

#### Scan Appliance    Photographs

Full Face / Full Smile + left / right / straight retracted in occlusion

